


**Tri-Tech
Medical Inc.**

*Manufacturer of
Medical Gas Pipeline Equipment*

Installation & Operating Instructions for Genesys ® 2 CC Series Manifolds






Features & Benefits


- Five year parts and one year labor limited warranty*
- Fully automatic changeover – no valves or levers to reset after each changeover
- Compatible with Tri-Tech Medical T-Net medical gas monitoring system saving you time and improving safety
- Field upgradeable design – kits allow unit to be changed from - i.e. cylinders to portable bulk or from standard flow to high flow or from lower delivery pressure to higher delivery pressure
- Field adjustable high/low line pressure alarms – using convenient buttons on circuit board.
- Circuit board triggers all required NFPA 99 alarms – simplifying wiring and reducing cost
- Unit includes hi/low line pressure transducer – eliminating need to purchase hi/low pressure switch, improving alarm accuracy and improving safety (manifold will automatically alarm if transducer goes bad)
- Easy to service layout/design
- Microprocessor based control panel incorporates LED's and illuminated text display readable even in poor lighting conditions
- Electronic monitoring of circuits, error, alarm or information messages displayed for ease of maintenance
- Accurate, long life pressure transducers for monitoring of line pressure and bank pressures
- Analog gauges also provided for use in event of power failure
- Pressures may be displayed in psi / kPa / BAR
- Built in DISS gas specific emergency feed ports
- Built in emergency reserve bank ports
- Input power 120 VAC, 50 to 60 Hz, (120-240 VAC 50 to 60 Hz on all models without heaters)
- Dual line pressure regulators on NFPA 99 models
- Gas specific header bar with integral check valves and cylinder pigtail assemblies
- Variety of header configurations available to meet the available space requirements of your installation
- Available in weatherproof cabinet for outdoor installation

* See Terms and Conditions, Document No. 99-0477, on our website at: www.tri-techmedical.com for complete details.



Contents

Features & Benefits -----	2	Appendix A	
General Instructions -----	4	Glossary of Terms -----	21
Introduction -----	4	Appendix B	
Installation-----	5	Technical Specifications -----	22
Mounting -----	5	Appendix C	
Plumbing -----	6	Model CC Piping Schematic-----	23
Electrical-----	8	Appendix D	
Remote Alarm Wiring-----	9	Operational Pressure Specifications --	24
Removing Transducer for PSM-----	10	Wiring Diagram – model CC-----	25
Installation -PSM-TXX-----	11	Wiring Diagram – model CC with heaters	
Installation – Pigtails & Cylinders ----	12	-----	26
Start Up and Checking Procedures ---	13	Wiring Diagram – Remote Audio/ Visual	
Start Up and Checking Procedures Cont.	14	Alarm TAV-1-----	27
Cylinder Replacement & Handling ---	15	Technical Assistance----	
Line Delivery Pressure Adjustment-----	15	Phone	800-253-8692
Programming Adjustments-----	16		or 440-937-6244
Table of Error Codes -----	17	Fax	440-937-5060
Table of Alarm Codes -----	17	E-mail	sales@tri-techmedical.com
Table of Information Codes-----	17		
Replacement Parts -----	18		
Trouble Shooting Guide-----	19		
Trouble Shooting Guide Cont.-----	20		
General Maintenance-----	20		



Introduction

Tri-Tech Medical manifolds are cleaned for use with oxygen. Each system is tested for changeover, triggering of alarms and leakage. Each unit is designed and prepared for the indicated gas service. Tri-Tech Medical manifolds are built in accordance with the National Fire Protection Association and Compressed Gas Association guidelines.

Warranty

All Tri-Tech Medical manifolds are warranted against defects in material and workmanship for the period of one year from date of purchase. All circuit boards are warranted against defects in material and workmanship for the period of three years from date of purchase.

General Instructions/Location & Shelter


Manifolds should be installed in accordance with guidelines stated by the National Fire Protection Association, the Compressed Gas Association, OSHA, and all applicable local codes. Central supply systems and cylinders should not be installed in a location where the temperature will exceed 125° F (51.6° C) or fall below -20° F (-29° C). A manifold placed in an open location should be protected against weather conditions. During winter, protect the manifold from ice and snow. In summer, shade the manifold and cylinders from continuous exposure to direct sunlight.

Leave all protective covers in place until their removal is required for installation. This precaution will keep moisture and debris from the piping interior.

Caution


Failure to follow the following instructions can result in personal injury or property damage:

- Never permit oil, grease, or other combustible materials to come in contact with cylinders, manifold, and connections. Oil and grease may react with explosive force when ignited while in contact with some gases – particularly Oxygen and Nitrous Oxide.
- Cylinder and master valves should always be opened very slowly. Heat of recompression may ignite combustible materials creating an explosive force.
- Pigtailed should never be kinked, twisted, or bent into a radius smaller than 3 inches. Mistreatment may cause the pigtail to burst.
- Do not apply heat. Oil and grease may react with explosive force when ignited while in contact with some gases – particularly Oxygen and Nitrous Oxide.
- Cylinders should always be secured with racks, chains, or straps. Unrestrained cylinders may fall over and damage or break off the cylinder valve which may propel the cylinder from its current position.
- Oxygen manifolds and cylinders should be grounded. Static discharges and lightning may ignite materials in an oxygen atmosphere, creating a fire or explosive force.
- Welding should not be performed near Nitrous Oxide piping. Excessive heat may cause the gas to dissociate, creating an explosive force.
- Remove all protective caps prior to assembly. The protective cap may ignite due to heat of recompression in an oxygen system.




Control Cabinet Installation

1. Determine and mark the vertical center line for installation of the manifold control cabinet.
2. If you wish to place cylinders under the manifold cabinet, measure from the floor to a point 79" in height above the finished floor of this vertical line. Using a level, mark a horizontal line at this point extending approximately 7" to the left and 7" to the right of center. This line indicates the location for the bottom of the Z mounting bracket. Mounting the Z bracket @ 79" aff* will result in the bottom of the manifold cabinet being 60" aff – allowing cylinders to be placed under the manifold cabinet. If you do not wish to place cylinders under the manifold cabinet, measure from the floor to a point 58 1/2" aff and follow the same steps above.
3. Mount the Z mounting bracket to the wall using fasteners suitable for the type of wall construction.
4. Temporarily hang the manifold cabinet on the Z bracket just installed, mating it with the Z bracket on the back of the manifold cabinet. Mark the locations of the two lower cabinet mounting holes on the wall.
5. Remove the manifold cabinet and install female portion of suitable fasteners for lower cabinet mounting holes.
6. Re-hang the manifold cabinet and install suitable fasteners in the lower cabinet mounting holes.




Header Installation


1. Attach the headers to the union on each side of the manifold control cabinet. Using a level, mark the placement of mounting brackets while keeping the header on a horizontal plane.
2. Remove the U – bolt assemblies from the header mounting brackets. Position the brackets so that the top of the bracket is aligned with the bottom of the headers and is centered between the cylinder connections. The end bracket should be placed as close to the last cylinder as possible to provide the most support and stability.



* add = above finished floor





1. Mark the mounting hole and install fasteners suitable for type of wall construction.
2. Fit the U – bolt over the header piping and tighten the two mounting nuts.




Plumbing – Model CC

1. The outlet of the manifold is located at the top center of the unit as shown here. The outlet connection is 1/2" NPT female. A source valve 1/2" M NPT x 3/4 tube extension has been provided with all NFPA 99 compliant models (part no. 48-0023) and should be installed in the manifold outlet. It may be necessary to install a 1/2" NPT 90° street elbow (part no.17-0300 sold separately) to meet the confines of the installation site.
2. Source valve part no. 48-0023
3. A 1/2" NPT brass pipe nipple has been provided (on CC models) to extend the intermediate relief piping outside the cabinet wall. The intermediate relief valve extension pipe nipple must be installed. It is bagged and shipped inside of the cabinet for protection during shipment.
4. It is also recommended that unions (part no. 17-0169) be installed between all relief valves and the relief (vent) pipeline system(s). There are two 1/2" NPT relief valve connections on the model CC manifold. See drawing at bottom of page 7.

Outlet 1/2" F NPT Line relief valve 1/2" F NPT





Plumbing – Model CC


5. A 1/2" NPT brass pipe nipple has been provided (on CC models) to extend the intermediate relief piping outside the cabinet wall.
6. The intermediate relief valve extension pipe nipple must be installed. It is bagged and shipped inside of the cabinet for protection during shipment.




Intermediate relief valve connection
1/2" NPT brass nipple provided
shown)

Recommended plumbing of relief and vent lines

(Items shown as dotted lines and unions are NOT supplied with CC manifold)






Electrical

1. Use one of the two 1/2" conduit knock-outs provided located nearest to the top left corner of the cabinet to route conduit to supply 120 VAC to the power supply. **Note: Separate conduit should be used for low voltage wires (use knock outs provided on the left side of the box).**

120-240 / 1 / 50 – 60 Hz may be used with all units EXCEPT those with heaters. Units with heaters may only be wired to 120/1/50-60. An additional transformer (sold separately part no. 35-3004) is required to connect units with heaters to 240 VAC



Conduit knockouts for
120 VAC

Conduit knockout for
low voltage alarm signals

2. Remove the power supply cover by loosening the screw located at the top of the cover and then sliding the power supply cover to the right until the screw is in the center of the tear-drop shaped cut out. Next, pull the cover forward until it clears the screw head and the fuse. Note: the bottom of the cover inserts into a slot in the back plate. Allow the cover to rest on the dual line regulator assembly plumbing just below the power supply.



3. Route wires of proper gauge (per local building code requirement) through the power supply conduit, thru the grommet on the power supply bracket and into the terminal strip.

Neutral

Load


Field Ground



4. Connect the 120 VAC facility **emergency power source** electrical wiring to the terminal strip provided on the front of the power supply mounting bracket (per photos right). (N = neutral, L = load, FG = field ground)



Note: The ground must be a solid earth ground with little or no resistance. A “noisy” earth ground may affect the digital display of the manifold.



Wiring – remote alarms

Caution: *Never connect or disconnect any electrical components with the power on. This may result in damage to components and is not covered under warranty*

1. Wires for remote alarms should be brought into the cabinet thru conduit or shielded cables (check local code requirements) thru the knockouts on the left side of the cabinet shown here. **Note: Separate conduit should be used for high voltage wires – never run low voltage wires in the same conduit as high voltage wires.**
2. If you are installing a model CC (cylinder x cylinder) cabinet there are three alarm signals required per NFPA 99, High Line Pressure, Low Line Pressure and Secondary in Use. The CC circuit board will trigger all three of these alarms (no hi/low pressure switch is required). The line pressure transducer must be installed outside of the cabinet – downstream of the source or main line valve with the cable being wired to the manifold circuit board to comply with NFPA 99. In this photo the line pressure transducer has been assembled into our PSM-TXX assembly and connected to the gauge port on the downstream (patient side) of the source valve and wired to the manifold circuit board.
3. The line pressure transducer may also be mounted inside the cabinet (as shown here). In this arrangement, a hi/low pressure switch (sold separately) will be required to meet the NFPA 99 recommendations. Note: the hi/low pressure switch would be wired directly to the master alarm panels – **not** to the manifold circuit board.
4. Remote alarm wires are connected to the circuit board at the terminal gate labeled MB-5. Signal wires and Common wires for Low Line Pressure, High Line Pressure and Secondary in Use should be connected to the terminals as indicated.
5. Note: all remote alarm terminals are normally closed when the gas pressure is in the normal range. The hi/low set points pre-programmed into the manifold circuit board logic chip are as per the charts on page 24. If desired, the high and low line pressure alarm set points may be adjusted by following the instructions on page 16.




Knock out for low voltage
remote alarm wiring



Line pressure transducer




MB-4						MB-5					
POWER SUPPLY INPUT						ALARM SIGNAL OUTPUTS					
ORANGE (+) 24 VDC	RED (+) 5 VDC	BLACK (-) COMMON	GREEN/EARTH GROUND	BLACK	BLACK	SOLENOID RIGHT BANK	SOLENOID LEFT BANK	RIGHT BANK/ SECONDARY LOW	LOW LINE PRESSURE	HIGH LINE PRESSURE	SECONDARY IN USE
				BLACK	BLACK				EMERGENCY RESERVE IN USE	EMERGENCY RESERVE LOW	
						NORMALLY CLOSED COMMON	NORMALLY CLOSED COMMON	NORMALLY CLOSED COMMON	NORMALLY CLOSED COMMON	NORMALLY CLOSED COMMON	NORMALLY CLOSED COMMON




Removing the Transducer and Demand Valve from inside the Manifold

Close the area zone valve. Warning: Do not close zone valve until all personnel have been advised of the intended service and all systems requiring medical gas are being supplied from portable or alternate supplies. Bleed the residual gas pressure from the manifold.


- 1 Carefully remove the Transducer from the Demand Valve. (*But do not discard the Transducer & Demand Valve.*) Clean off any debris from the threads of the transducer nut and nipple.




- 2 Carefully remove the Demand Valve from the Outlet Block and remove any loose debris from the valve. Also carefully remove any loose debris from the threads of the Outlet Block.




- 3 Once the debris is removed from Outlet Block, replace the Demand Valve with a 1/4" NPT Hex Plug that comes with the kit. Thread on and tighten the plug onto the Outlet Block.






Installing the PSM-T Valve Kit outside the Manifold

1




Remove the hex plug from the Source Valve side port. Carefully remove any debris from the threads of the Source Valve.

2




Select one of the two adapters that come with the PSM-T Valve Kit. After selecting the proper adapter, thread on and tighten the adapter onto the Source Valve side port.

3




With the proper adapter selected. Re-wrap the threads of the Demand Valve with oxygen clean wrapping tape. Thread on and tighten onto the Source Valve. *(This is same Demand Valve from inside the manifold cabinet).*

4




With the Demand Valve tighten, thread on and tighten the PSM-T on the Demand Valve.

5




With the PSM-T is tighten down, re-insert the Transducer through the manifold. *(This is the same Transducer from inside the manifold cabinet).* Thread on and tighten the Transducer to the bottom of the PSM-T valve Kit.

6




The PSM-T Valve kit is completed



Installing pigtails & attaching cylinders – CC models

1. The check valve outlet fittings on the manifold header bars are CGA (Compressed Gas Association) gas specific threads. Each of these fittings has an integral check valve. Make sure the 3 Digit CGA number stamped on the outer perimeter of these fittings matches the CGA number stamped on the mating CGA fittings on the pigtails.
2. Connect the pigtails to the check valve outlets on the manifold headers.
3. Check the master valves to be certain they are open (turn counter-clockwise to open). (Note: the master valve should always be left open. It is to be used only in the event of an emergency).
4. SLOWLY open all cylinder valves (turn counter-clockwise to open). Check all cylinder and pigtail connections for leaks using an oxygen safe leak test solution (any bubbles forming around connections indicate leakage).
5. Note: The manifold has been tested for leaks at the factory, but the installer MUST check for leaks at all connections made during installation.


Check valve
outlet fittings




Line transducer
mounted on patient
side of source valve
using PSM-TXX
xxx assembly.

Master valves

Attach pigtails to
header check valve
outlets using 1-1 /8”
open end wrench







Start up & checking procedures


The manifold is pre-programmed (per page 24) and tested before it is shipped. You may, however, wish to modify some of the programming (see page 16). The unit has been designed to allow some programming to be simply and safely altered in the field.

1. Start with all cylinders turned off and with zero pressure supplied to the manifold cabinet. To conduct the initial start-up testing of the manifold, it is simpler and faster if the circuit board is switched from the standard 'cycling view' (this is the mode in which it is shipped) to the 'global view'. To do this, the cover must be removed from the circuit board and the top No.1 dip switch must be switched from the RUN mode to the PROGRAM mode. (**Note:** local alarms are **not** displayed when in the global view mode).
2. Turn on the 120 VAC to the unit. For model CC the display should illuminate showing all zeros for the Line Pressure (LIN), Left Bank Pressure (LEFT) and Right Bank Pressure (RIGHT). The intermediate (INT) and emergency reserve (EMER) displays should also display zeros. The INT and EMER displays will always display zeros on the global view screen of CC models and will not display when the mode is switched to the cycling view. (CC models). The INT and EMER displays are only activated on LL models. Both the left & right bank Red (Depleted) LED's should be illuminated. Both the left and right bank Green (In Use) and Yellow LED's (Ready) should be extinguished.
3. SLOWLY open one-cylinder valve on the left bank. The left bank pressure gauge (inside the cabinet) and the text display (on the outside of the cabinet) should show the full pressure of the cylinder. The Red (Depleted) LED for the left bank should have extinguished leaving only the Green (In Use) LED illuminated.
4. SLOWLY open one-cylinder valve on the right bank. The right bank pressure gauge (inside the cabinet) and the digital display (on the outside of the cabinet) should show the full pressure of the cylinder. The Red (Depleted) LED for the right bank should have been extinguished and the Yellow (Ready) LED should have illuminated.

Program control buttons


Dip switch #1 RUN/PROGRAM






Start up & checking procedures (cont'd)


5. Turn off all open left bank cylinder valves. Create a slight flow of gas in the delivery pipeline system. DISS demand valves have been provided on the line regulators. Mating DISS fittings may be used to create a flow of gas within the manifold cabinet. The left bank pressure text display and pressure gauge should fall and the control automatically switches over to the right bank. Delivery pressure remains constant. The left bank Red (Empty) LED will illuminate. The secondary supply in use alarm should activate on the master alarm(s).




6. SLOWLY reopen the cylinders on the left bank. The left bank pressure gauge and digital display should return to full pressure. The left bank yellow (Ready) LED will illuminate simultaneously the left bank red (Empty) will extinguish. All remote secondary supply in use alarms will be canceled.




7. Repeat steps 5 & 6 to simulate an empty right bank.
8. After you are satisfied that the manifold is functioning properly and that all master alarm signals are being triggered properly, the manifold circuit board should be returned to the 'cycling view' mode. This is achieved by moving the dip switch from the PROGRAM to the RUN position.






Cylinder replacement & handling


1. Close all cylinder valves on the depleted bank.
2. SLOWLY loosen and remove the pigtail connection from the depleted cylinders.
3. Remove depleted cylinders and replace protective caps.
4. Place and secure full cylinders into position using chains, belts or cylinder stands.
5. Remove protective cylinder caps from full replacement cylinders. With the valve outlet pointed away from all people in the area, slowly open each cylinder valve slightly for a split second to blow out any dirt or contaminants that may have become lodged into the cylinder valve.
6. Connect pigtails to cylinder valves and tighten with wrench.
7. SLOWLY turn each cylinder valve until each cylinder is fully on.
8. Leak test the connections using an oxygen approved leak test solution.
9. Observe the following conditions: The red (Empty) LED is extinguished and the yellow (Ready) LED is illuminated and the secondary in use alarm is cancelled.
10. The manifold supply bank is now replenished and automatically placed in “reserve”.



Line delivery pressure adjustment

1. Leave the manifold in full operational status.
2. Create a flow condition in the delivery piping system. DISS demand valves have been provided on the line regulators for this purpose. Mating DISS fittings may be used to create a flow of gas within the manifold cabinet.
3. Open the manifold cabinet door and locate the line pressure regulators. Ball valves on the inlet and outlet sides of each regulator determine which regulator is “on line” and which is “off line”. Note: when the ball valve handle is perpendicular to the pipeline, the ball valve is closed. One-line regulator should be valve closed and the other valve open. The line regulators on-line and off-line should be alternated every other month to ensure each is exercised.
4. Turn the T – bar handle clockwise to increase pressure or counter-clockwise to decrease pressure. It may be necessary to use a 3/4” open-end wrench, loosen the locknut on the adjusting screw (*on High Flow Models Only*).
5. After adjustment, retighten the locknut on the adjusting screw and close the cabinet door.





Programming Adjustments

The manifold is pre-programmed and tested before it is shipped. You may, however, wish to modify some of the programming. The unit has been designed to allow some programming to be simply and safely altered in the field. The aspects of the program which may be altered include:



- ✓ Low line pressure set point
- ✓ High line pressure set point
- ✓ Units of measure (psi, bar or kPa)
- ✓ Calibrate line sensor
- ✓ Display scroll time (seconds)
- ✓ Application logic (this should only be changed if the unit is being converted to a new application – in conference with a factory technician)


1. To make any programming change, the cover must be removed from the circuit board and the top #1 dip switch must be switched from the RUN mode to the PROGRAM mode.
2. Once the #1 dip switch has been changed from the RUN mode to the PROGRAM mode, the display should look like this
3. While viewing the display on the front of the cabinet door, use the NEXT or PREVIOUS buttons to toggle thru the menu choices. When you find the item you wish to reprogram, use the UP or DOWN buttons to display the new setting desired and then use the ENTER button to save changes to the programming. **Note – if the ENTER button is not pressed for each and every change, that programming change will not be saved.** The display will display the word SAVED when you have successfully saved a programming change.
4. Return the dip switch to the RUN position and replace the circuit board cover when you are finished making changes to the programming.

Displays showing common pre-programmed alarm set points

Program control buttons

Dip switch #1 RUN/PROGRAM





Error Codes, Alarm Codes & Information Codes


Code Displayed	Message Displayed	Explanation
Error 01	Left bank sensor out of range	This condition is activated when the left sensor's raw readings are at either extreme. Can be caused by a disconnected, wired incorrectly, bad, or over-pressurized sensor.
Error 02	Right bank sensor out of range	This condition is activated when the right sensor's raw readings are at either extreme. Can be caused by a disconnected, wired incorrectly, bad, or over-pressurized sensor.
Error 05	Left bank pressure high	In CC models this message is displayed whenever the left inlet bank pressure exceeds 3,000 psi. In LL models this message is displayed whenever the left inlet bank pressure exceeds 500 psi.
Error 06	Right bank pressure high	In CC models this message would be displayed if the right inlet bank pressure exceeds 3,000 psi. In LL models this message would be displayed if either inlet bank pressure exceeds 500 psi.
Error 09	Line sensor noise detected	Message is displayed if the circuit board detects noise in the signal from its digital sensor. Noise is detected if the protocol checksums do not match.
Error 10	Line sensor failed to respond	Message is displayed if a line sensor is not responding.
Error 11	Line sensor is disconnected	Message is displayed whenever a line sensor is disconnected.
Error 12	Secondary supply leak detected	Message is displayed when a leak is detected in the secondary bank. (CC models only).
Error 14	Gas type mismatch	This error code is triggered by a mismatch in gas type between the line sensor and user selected gas type in setup of the manifold circuit board.
Alarm 01	Line pressure low	Message is displayed and low line pressure relay activated whenever the line pressure is below the programmed low line pressure alarm set point.
Alarm 02	Line pressure high	Message is displayed and high line relay activated whenever the line pressure is above the programmed high line pressure alarm set point (high line pressure alarm is triggered).
Alarm 03	Secondary supply in use	Message is displayed and secondary in use relay activated when the manifold has switched over to the secondary bank. Clears when tank is replaced.
	XX% Remains	(CC models only - except for N2O & CO2 services). It calculates the percent of gas remaining in the service bank.

Definitions and clarifications


Alarm Code – Alarm conditions per NFPA 99C and Z7396.1 guide lines.

Error Code – Messages that provide diagnostic information to assist in resolving system problems.

Information Codes – Messages that provide information regarding the operation




Item	Part Number	Description
Primary Regulator & Repair Kit	68-0003R	Primary regulator
	68-0003RK	Primary regulator repair/rebuild kit
Line Regulators & Repair Kits	68-0004R	Line regulator standard flow 5 - 125 psi
	68-0004RK	Standard flow line regulator repair/rebuild kit
	68-0002R	Line regulator high flow 5 - 125 psi
	68-0002RK	High flow line regulator repair/rebuild kit 5 - 125 psi
	68-0001R	Line regulator high flow 10-200 psi
	68-0001RK	High flow line regulator repair/rebuild kit 10 - 250 psi
Circuit Board	35-1003R	CC series PLC board with single text display
Power Supply	35-2013R	Power supply
Transducers/Sensors	14-3001R	0 - 2,500 psi transducer with 3' cable for left or right banks
	14-3024	0 - 250 psi transducer with 1.5' cable N2
	14-3025	0 - 100 psi transducer with 1.5' cable Oxygen
	14-3026	0 - 100 psi transducer with 1.5' cable Med Air
	14-3027	0 - 100 psi transducer with 1.5' cable N2O
	14-3028	0 - 100 psi transducer with 1.5' cable CO2
Solenoid Valves	48-1007R	Solenoid valve for CC series
Check Valve	17-4003R	Intermediate check valve 1/2" M NPT x 1/2" OD tube
Tubes & Compression Fittings	17-4012	Compression Sleeve 1/2" OD tube - glass filled Teflon
	17-4005	Compression Nut for 17-4012
	Q1100-1	1/2" OD copper tube x 7"
	17-4013	Compression Sleeve 3/8" OD tube - glass filled Teflon
	17-4024	Compression Nut for 17-4013
Gauges	14-1018	0 - 4,000 psi 1 1/2" x 1/8" M NPT center back gauge
	14-1016	0 - 400 psi 2" x 1/4" M NPT bottom port gauge
	14-1017	0 - 400 psi 1 1/2" x 1/8" M NPT center back gauge
	14-1009	0 - 300 psi 1 1/2" x 1/8" M NPT center back gauge
	14-1008	0 - 100 psi 1 1/2" x 1/8" M NPT center back gauge
Relief Valves	RV-22-075	75 psi x 1/2" M NPT inlet with pipe away adaptor
	RV-22-150	150 psi x 1/2" M NPT inlet with pipe away adaptor
	RV-22-250	250 psi x 1/2" M NPT inlet with pipe away adaptor
	RV-11-400	400 psi x 1/4" M NPT inlet with pipe away adaptor
Pigtails	20-1001	36" (pre-bend) single loop rigid copper O2 (w/o check) – CGA 540
	20-1002	36" single loop rigid copper N2O – CGA 326
	20-0002	24" Flexible stainless braided N2O - CGA 326 with captured fittings
	20-0003	24" Flexible stainless braided CO2 – CGA 320 with captured fittings
	20-0004	24" Flexible stainless braided AIR – CGA 346 with captured fittings
	20-0005	24" Flexible stainless braided N2 – CGA 580 with captured fittings
Union	17-0169	Union 3-piece 1/2" F NPT x 1/2" F NPT & (1) 5" lg 1/2" M NPT Nipple
Master Valve	17-0256	Master Valve 1/2" F NPT x 1/2" F NPT
Master Valve Repair Kit	GMV-1001RK	Master valve rebuild kit
Heater element	35-2001	Ceramic Heater
Accessory	TAV-1	Manifold Remote Audio/Visual Alarm kit



Note: trouble-shooting and repairs should be done by qualified personnel only.

Component	Symptom	Probable Causes	Remedy or Check
Circuit board	No indicator lights on front panel illuminate when power is connected.	Wiring connection. Power input out of range Bad circuit board	Check wiring connections. Check electrical power supply output voltages. Replace power supply or circuit board.
Circuit board	Red indicator lights are on but both banks are full	Master valve or cylinder valves on bank are closed. Pigtails are installed with check valves in wrong direction Bank pressure is not sufficient for logic board to place it in "In Use" or "Ready" status (see appendix E minimum inlet pressures page 27)	Open valves SLOWLY Close cylinders and re-install pigtails in proper flow direction Replace cylinders with full cylinders. Or, if using portable bulk vessels, open pressure building valve on vessel or replace portable bulk vessel with higher delivery pressure portable bulk vessel
Circuit board	Error code(s) being displayed	Loose or disconnected or broken wire, mis-connected wire, a bad transducer, a calibration problem or an over-pressure situation.	Check wires for good and correct location connection to circuit board (See ERROR CODES page 17) If all wires are connected properly and located properly – it may be necessary to replace a transducer.
Entire manifold	Loss of cylinder contents	Leakage in manifold cabinet, headers or pigtails. Leakage thru manifold solenoid vent / relief Leakage around regulator bonnet Regulator with bad seat Leaking gauge Regulator set too high Overpressure due to failed regulator seat Regulator freeze-up (N2O or CO2) / heater failure	Locate leak using oxygen compatible leak test solution, tighten, reseal or replace leaking fitting(s) or pigtails. Replace solenoid valve Tighten regulator bonnet Rebuild or replace regulator Replace gauge Set delivery pressure per specifications (see page 24) Rebuild or replace regulator Flow rate exceeds manifold design capacity or cylinder withdrawal capacity – add cylinders. Repair heater or add heater and consider adding additional cylinders




Note: trouble-shooting and repairs should be done by qualified personnel only.

Component	Symptom	Probable Causes	Remedy or Check
Loss of cylinder contents	Both banks feeding	Leaking header/pigtail connection	Tighten fitting or re-tape with Oxygen safe Teflon tape (if NPT fitting) and tighten.
		Leaking intermediate check valve	Replace check valve
		Leaking solenoid valve	Replace solenoid valve
		Primary regulator(s) out of adjustment	Set delivery pressure to specifications per chart on page 24
System	Changeover occurs, secondary in use alarm is triggered and then clears placing empty or near empty bank into reserve	Manifold is unable to support required flow (N2O or CO2 units)	Increase manifold flow capacity (call factory for assistance) Increase bank size

General Maintenance

	Control Cabinet	Headers & Pigtail
Daily	Record line and bank pressures	Observe nitrous oxide and carbon dioxide systems for cylinder frosting or surface condensation. Should excessive condensation or frosting occur it may be necessary to increase manifold capacity.
Monthly	Check regulators, compression fittings and valves for external leakage. Check valves for closure ability. Alternate line regulator in use (if dual). Check regulators for stable delivery pressure.	Inspect valves for proper closure. Check pigtails for cleanliness, flexibility, wear, leakage, kinked, pinched or twisted and thread damage. Replace damaged pigtails immediately. Inspect header check valve outlets for closure ability.
Annually	Check relief valve pressures Check regulator seats.	Check pigtails for cleanliness, flexibility, wear, leakage, kinked, pinched or twisted and thread damage. Replace damaged or worn pigtails immediately.
	Tighten regulator bonnets approximately 1 to 2 degrees (out of 360).	
Every 4 years		Replace all pigtails
Every 10 years	Rebuild or replace: all regulators, 17-4003 intermediate check valves, 48-1007 solenoid valve, and DISS demand check valves.	Rebuild or replace master valves. Replace all check valve outlets.



Appendix A

Glossary of Terms

AC *Alternating Current*

An electric current that reverses direction or polarity at regular intervals.

Alarm Code

Alarm conditions per NFPA 99C guidelines.

BAR *Bar*

A measurement of the force in a compressed gas system.

1 BAR = 14.7 psi (1 atmosphere)

Check Valve

A valve which operates mechanically and automatically to stop a reverse flow of gas

DC *Direct Current*

An electric current that flows in one direction. The current can be steady or pulse.

Economizer Circuit

A mechanical piping circuit which allows built up reserve gas to be used in low volume rather than allowing the gas to vent to atmosphere.

Error Code

Messages that provide diagnostic information to assist in resolving system problems.

Information Code

Messages that provide information regarding the operation of the system.

kPa *Kilopascals*

A measurement of the force in a compressed gas system.

1 kPa = .14 psi

LED *Light Emitting Diode*

A semiconductor diode that converts applied voltage to light.

NFPA *National Fire Protection Association*

The National Fire Protection Association is an association engaged in standards development.

NO *Normally Open*

An electrical circuit in which the switch is normally open. No current flows through the circuit in normal operation. Only when the switch is closed is the flow of current started.

A normally open solenoid valve is one designed so that it is open when there is no power to the solenoid and closed when it is energized.

NC *Normally Closed*

An electrical circuit in which the switch is normally closed. Current flows through the circuit in normal operation. Only when the switch is opened is the flow of current stopped.

psi *Pounds per Square Inch*

A measurement of the force in a compressed gas system.

1 psi = 6.9 kPa

Solenoid Valve

A valve that is opened or closed electromagnetically.

Transducer


A device that converts pressure into an electrical signal.

Transient Signal

An intermittent and brief signal that quickly corrects and returns the alarm to a normal operating mode before monitoring personnel can silence the alarm

V *Voltage*

Voltage is electrical pressure or force. One volt is equal to the difference of electrical potential between two points on a conducting wire carrying a constant current of one ampere when the power dissipated between the points is one watt.



Appendix B – Technical Specifications

Maximum Inlet Pressure: Models CC: 3,000 psi

Operating Ambient Temperature range:

Model CCU: -20 F (-29 C) to 125 F (54.4 C) all gas services

Storage Temperature: -20C (-4F) to + 185F (85 C)

AC Input: 120 volts AC - 50-60 Hz

Input Fuse: 5 Amp input AC line fuse protects the input wiring to power supply

Power Consumption: 45W (0.4 amps using 120 VAC) maximum without heaters
245W (2.1 amps using 120 VAC) maximum with heaters

Pressure Measurement Accuracy:

0-100 psi transducer +/-1% of full scale - Line Pressure
Oxygen, Nitrous Oxide, Medical Air, Carbon Dioxide

0-250 psi transducer +/-1% of full scale - Line Pressure - Nitrogen

0 – 2,500 psi transducer +/-2% of full scale – Bank Pressures

Solenoid Valve: 24VDC – Normally Open (Valve opens when de-energized)

Dimensions

Control Cabinet: Dimensions excluding inlet & outlet fittings
15 3/8" W x 25" H x 9 1/4" D



Dimensions including inlet & outlet fittings
17 1/4" W x 27" H x 9 1/4" D


Line Pressure Transducers: Housing dimensions: 1.25" Diameter x 3.75" Length
including inlet fittings

Weatherproof Control Cabinet

Dimensions excluding inlet & outlet fittings
17 1/4" W (cabinet) 18 3/4" W (door) x 26 3/4" H x 11" D

Dimensions including inlet & outlet fittings
20 1/4" W x 29" H x 11" D





Appendix D – Operational Pressure Specifications

(all pressures shown in psi)

Minimum inlet pressure requirements for CC Manifold


Manifold delivery pressure	Minimum inlet pressure Left bank	Minimum inlet pressure Right bank
50	300	250
80	300	250
170	300	250

Alarm Pressure Settings for CC Models



Normal Delivery Press	Line Relief Setting	High Line Press Set Point	Low Line Press Set Point	Secondary in Use Set Point
50	75	60	40	300 left bank / 250 right bank
80	150	96	64	300 left bank / 250 right bank
170	250	200	140	300 left bank / 250 right bank


Primary Regulator Settings (CC Models) Note: All settings done with full cylinder pressure and with slight gas flow thru the manifold. DISS demand valves have been provided on the line regulators and may be used to create a slight flow. Primary regulator outlet pressure will vary with varying inlet pressures. (The outlet pressure will rise up to 30 psi above set point as the cylinder pressure decreases).

Normal Delivery Pressure (factory delivery pressure)	Left Primary Regulator Set Point	Right Primary Regulator Set Point
50	260	220
80	260	220
170	260	220



[illegible]





**Tri-Tech
Medical Inc.**

*Manufacturer of Medical
Gas Pipeline Equipment*

This page is left Intentionally Blank



**Tri-Tech
Medical Inc.**

Tri-Tech Medical Inc., 35401 Avon Commerce Pkwy., Avon, Ohio 44011

Tel. 1-800-253-8692 -or- 440-937-6244 Fax. 440-937-5060

Web site www.tri-techmedical.com E mail address: sales@tri-techmedical.com